



FRIENDSHIP MEDICAL



FRIENDSHIP MEDICAL

西安富德医疗电子有限公司
Xi'an Friendship Medical Electronics Co., Ltd.

Electroneurodiagnostic
Electrodes and Accessories

- Friendship Medical is engaged to providing the customers with high-quality and low-cost Disposable Subdermal Needle Electrodes, Disposable Pre-Gelled Ag/AgCl Surface Electrodes, Disposable Nerve Stimulator Probes and other electroneurodiagnostic accessories.
- Meeting your requirements is the cornerstone of our customer service policy.
- Our products have got the certification of American FDA, European CE, Canadian CMDCAS, MSA in Japan and KFDA in Korea!
- We can produce various kinds of complicated and cutting-edge medical products under OEM mode!
- We are an OEM company with large production capacities.
- Our company has 5 categories, 18 series and 268 models of electrodes. The annual output of different kinds of electrodes is 20,000,000 sets.

西安富德医疗电子有限公司
Xi'an Friendship Medical Electronics Co., Ltd.

Version H / August 18, 2013



Headquarters / 中国总部
Add / 地址: No.9 Gao Xin 1st Road, Hi-Tech Development Zone, Xi'an, Shaanxi, China, 710075
Tel / 电话: 86 29 88225200
Fax / 传真: 86 29 88236285
URL / 网址: www.xafdec.com
Email / 邮箱: georgezhai2616@163.com

美国分公司 / USA Branches
Los Angeles: Max Neuro Supply Inc.
17565 E. Walnut Dr.
South, City of Industry, CA 91748
Tel: 626-757-5513
Fax: 626-270-4050
email: edward@maxneurosupplyinc.com

美国分公司 / USA Branches
Maryland: ZYTEK Neuro Inc.
8921 Liberty Lane, Potomac,
MD 20854 USA
Tel: 240-838-4498
Fax: 301-251-8718
www.zytekneuro.com
email: ssliu87@hotmail.com

加拿大分公司 / Canada Branches
Winnipeg: CareTrax Neuro Inc.
3 St. Anne's Road, Winnipeg,
Manitoba R2M 2X9, Canada
Phone: 204-926-8546
Fax: 204-926-8506
email: office@caretraxneuro.com

欧洲分公司 / Europe Branches
Denmark: Friendship Europe ApS
GI.Landevej 20
4000 Roskilde, DENMARK
Tel: +45 20164411
email: erik.jensen@friendshipeurope.com

Disposable Subdermal Needle Electrode



Disposable Subdermal Needle Electrodes, Single, Multicolor



Disposable Subdermal Needle Electrodes, Twisted Pair, Multicolor



Disposable Subdermal Corkscrew Needle Electrodes (Single, Twisted, Multicolor)



Disposable Subdermal Needle Electrodes Short Lead Length in 10cm Needle with Short Lead Length bent in 36°



Disposable Subdermal Needle Electrodes (Parallel Pair, Multicolor)



Disposable Subdermal Needle Electrodes (Five / Six Parallel)



Medical Extension Lead Wire for Disposable Subdermal Needle Electrodes (with Short Lead Length)



Disposable EMG Concentric Needle Electrodes



Disposable EMG Monopolar Needle Electrodes



5-Pin-DIN Plug Cable for Monopolar Needle Electrodes



Cable for Concentric Needle Electrodes



Disposable Botox Needle Electrodes



Connection Method of Disposable EMG Monopolar Needle Electrodes and Single Cable



Disposable EMG Sensory Needle




Connection Method of Disposable EMG Monopolar Needle Electrodes and the 5-Pin-DIN Plug Cable




High-Field MRI Intraoperative Electroneurodiagnostic Specific Probes and Electrodes (for Eliminate MRI Artifacts and Image Distortion)




Disposable Pre-Gelled Ag/AgCl Surface Electrode




Disposable Pre-Gelled Ag/AgCl Surface Electrodes (Single)




Disposable Pre-Gelled Ag/AgCl Surface Electrodes (Twisted-pair)



Disposable Central Ag/AgCl Surface Electrodes (Use with Mini-Crocodile Clip)




Shielded Electrode Cable




Shielded Electrode Cable




Crocodile Clip Cable



Shielded Electrode Cable



Disposable Central Ag/AgCl Surface Electrodes (with Lead Wire)



Disposable Pre-Gelled Ag/AgCl Surface Electrodes (Connector: 0.7mm, Single)



Disposable Central Ag/AgCl Surface Electrodes (with Lead Wire)




Disposable Pre-Gelled Ag/AgCl Surface Electrodes (Round, Single)


Disposable Subdermal Needle Electrode




Disposable Concentric Direct Nerve Stimulator Probe




Disposable Ball Tip Direct Nerve Stimulator Probe (90mm-300mm)




Disposable Flush Tip Direct Nerve Stimulator Probe




Disposable Monopolar Direct Nerve Stimulator Probe



Disposable Hook Nerve Stimulator Probe (Triple Hook 100° / 180°)



Disposable Hook Nerve Stimulator Probe (Double Hook 100° / 180°)



Disposable Double Ball Tip Direct Nerve Stimulator Probe



Disposable Bipolar Direct Nerve Stimulator Probe



Disposable Ball Tip Direct Nerve Stimulator Probe (90mm-150mm)



Various Double Hook Nerve Probes
Double Hook 100° / Double Hook 180°
Peacock Hook 90° / Double Ball Tip

Medical EEG Cup Electrode



Medical EEG Gold Cup Electrodes
(Lead Wire Contains Two Kinds of
Wires: Long-Life Teflon Wire
and Silicone Wire)



Medical EEG Ag Cup Electrodes
(Lead Wire Contains Two Kinds of
Wires: Long-Life Teflon Wire
and Silicone Wire)



Medical EEG Ag/AgCl Cup Electrodes
(Lead Wire Contains Two Kinds of
Wires: Long-Life Teflon Wire
and Silicone Wire)



Din Jack



Extension Lead Wires (4 Colors)



Disposable Digital
Ring Electrodes



Harness for 18 Channels
of Disposable
Subdermal Needle Electrodes

